

Whole Mount smFISH with HCR Protocol

Optimized for *Aurelia* jellyfish, modified from Choi *et al.* (2014) ACS Nano, 8(5): 4284–4294

Last Updated 7/25/17

Document created by David Gold (www.DavidAdlerGold.com)

Antisense probe(s)	Sense probe(s)

Fixation:

- 1. Anesthetize animals in 7.3% MgCl₂
- 2. Fix animals in 4% formaldehyde for 1 h @ room temperature (RT).

Acclimate from 7.3% MgCl₂ to MeOH using 5 minute washes:

- 3. 25% MeOH
- 4. 50% MeOH
- 5. 75% MeOH
- 6. 100% MeOH
- 7. Store animals at -10°C

Rehydration:

5 minute washes at room temperature:

- 8. 60% MeOH 40% PBSTr
- 9. 30% MeOH 70% PBSTr
- 10. 100% PBSTr
- 11. Wash 3x in PBSTr over a 2-hour period at room temperature (RT): (1) (2) (3)

Two-stage multiplexed in situ hybridization using DNA HCR Detection stage:

- 12. Pre-hybridize with 350 μ L of probe hybridization buffer for 30 min at 45°C.
- 13. Prepare probe solution by adding 1 pmol of each probe (1 μ L of 1 μ M stock per probe) to 500 μ L of probe hybridization buffer at 45°C.
- 14. Remove the pre-hybridization solution and add the probe solution.
- 15. Incubate the embryos overnight (12–16 h) at 45°C.

Remove excess probes by washing at 45 °C with 500 µL:

- *Wash solutions should be pre-heated to 45 °C before use.*

- 16. 75% of probe wash buffer / 25% 5x SSCT for 15 min
- 17. 50% of probe wash buffer / 50% 5x SSCT for 15 min
- 18. 25% of probe wash buffer / 75% 5x SSCT for 15 min
- 19. 100% 5x SSCT for 15 min
- 20. 100% 5x SSCT for 30 min.

Amplification stage

- 21. Pre-amplify embryos with 350 µL of amplification buffer for 30 min at RT.

- 22. Prepare 30 pmol of each fluorescently labeled hairpin by snap cooling in 10 µL of 5x SSC buffer (heat at 95 °C for 90 seconds and cool to room temperature on the benchtop for 30 min).

- 23. Prepare hairpin solution by adding all snap-cooled hairpins to 500 µL of amplification buffer at room temperature.

- 24. Remove the pre-amplification solution and add the hairpin solution.

- 25. Incubate the embryos overnight (12–16 h) at room temperature.

Remove excess hairpins by washing with 500 µL of 5x SSCT at RT:

- 26. 2 x 5 min: (1) (2)
- 27. 2 x 30 min: (1) (2)
- 28. 1 x 5 min

RECIPIES (Make buffers fresh every 2-3 weeks)

PBSTr (50mL)

150 µl Triton X-100 in 50mL PBS

Probe hybridization buffer

For 40 mL of solution

50% formamide

20 mL formamide

5x sodium chloride sodium citrate (SSC)

10 mL of 20x SSC

9 mM citric acid (pH 6.0)

360 µL 1 M citric acid, pH 6.0

0.1% Tween 20

400 µL of 10% Tween 20

50 µg/mL heparin

200 µL of 10 mg/mL heparin

1x Denhardt's solution

800 µL of 50x Denhardt's solution

10% dextran sulfate

8 mL of 50% dextran sulfate

Fill up to 40 mL with ultrapure H₂O

Probe wash buffer

For 40 mL of solution

50% formamide

20 mL formamide

5x sodium chloride sodium citrate (SSC)

10 mL of 20x SSC

9 mM citric acid (pH 6.0)

360 µL 1 M citric acid, pH 6.0

0.1% Tween 20

400 µL of 10% Tween 20

50 µg/mL heparin

200 µL of 10 mg/mL heparin

Fill up to 40 mL with ultrapure H₂O

Amplification buffer

For 40 mL of solution

5x sodium chloride sodium citrate (SSC)

10 mL of 20x SSC

0.1% Tween 20

400 µL of 10% Tween 20

10% dextran sulfate

8 mL of 50% dextran sulfate

Fill up to 40 mL with ultrapure H₂O

5x SSCT

For 40 mL of solution

5x sodium chloride sodium citrate (SSC)

10 mL of 20x SSC

0.1% Tween 20

400 µL of 10% Tween 20

Fill up to 40 mL with ultrapure H₂O

50% dextran sulfate

For 40 mL of solution

50% dextran sulfate

20 g of dextran sulfate powder

Fill up to 40 mL with ultrapure H₂O